FANTESTIC[®] Z CORE DC

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FANTESTIC® Z CORE DC

CE 0482

Hersteller

R-dental Dentalerzeugnisse GmbH • Winterhuder Weg 88 • 22085 Hamburg
 Telefon 0 40 - 30707073-0 • Fax 040 - 30707073-73 • www.r-dental.com

FANTESTIC® Z CORE DC (KwikkMix®, 5 ml)**Bestell-Nr.: FCS4065 (Farbe A2 Universal)****Bestell-Nr.: FCS4067 (Farbe weiß-opak)****FANTESTIC® Z CORE DC (AutoMix, 25 ml)****Bestell-Nr.: FCK4065 (Farbe A2 Universal)****Bestell-Nr.: FCK4067 (Farbe weiß-opak)****Beschreibung und Indikation**

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FANTESTIC® Z CORE DC is a self-curing, fluoride filler containing, radiopaque core build-up composite (hybrid composite) with light-curing option (dual-cure, DC).

FANTESTIC® Z CORE DC is indicated for adhesive core build-ups of vital and non-vital teeth, fillings by insufficient tooth structure and for adhesive cementing of glass fibre reinforced composites (e. g. pre-silanized glass-fiber posts, like FANTESTIC® POST).

FANTESTIC® Z CORE DC contains nano-zirconium dioxide and nano-calciumfluoride fillers. These lead to a high compression strength and a strong radiopacity. FANTESTIC® Z CORE DC can be cut like dentin. Because of the nano-zirconium and nano-calciumfluoride fillers an excellent standability, consistency and wetting of tooth structure ("wettability") is achieved.

On teeth structure FANTESTIC® Z CORE DC is used with dual-curing or self-etching, dual-curing adhesive systems (see below 3 b. FANTESTIC® FLOWSIVE or FLOWSIVE 2 with ACTIVATOR or 4. FANTESTIC® FLOWSIVE SE).

Root canal posts made of fiber reinforced composites have to be silanized with FANTESTIC® CerBond before luting.

Root canal posts made of metal, all-ceramics or zirconium dioxide could be primed with the universal primer ZM-Bond®. Pay attention to the manufacturer's instructions.

Characteristics and composition

FANTESTIC® Z CORE DC is based on multifunctional methacrylates and inorganic filler particles like barium glass, nano-zirconium dioxide and nano-calciumfluoride. The filler content is about 70 % by weight.

FANTESTIC® Z CORE DC is delivered in commercial 5 ml-double-chamber syringes (KwikMix®) and 25 ml-cartridges (AutoMix).

Working time and hardening time

Working time and hardening time of FANTESTIC® Z CORE DC depends on temperature. Immediately after automatically mixing:

Extraoral application at about 23 °C ± 1 °C (73 °F ± 2 °F) until insertion* intraoral at about 37 °C ± 1 °C (99 °F ± 2 °F):

Working time	1 min 30 s ± 10 s
Hardening time (without working time)	3 min 30 s ± 10 s

*with FANTESTIC® FLOWSIVE / FLOWSIVE 2 with ACTIVATOR or with FANTESTIC® FLOWSIVE SE

FANTESTIC® Z CORE DC meets the requirements of ISO 4049.

Recommended procedures

1. Isolation

A dry area and use of rubber dam (Order-No.: KDL31XX, XX = depends on strength) or latex-free dam (Order-No.: KDL3500) is recommended.

2 a. Preparation

Prepare in accordance with the clinical rules, disinfect and dry and in case the pulp should be protected with a Ca(OH)₂ cement (for example Cal de Luxe®, Order-No.: CDL2503).

2 b. Preparation of the root canal

Prepare the root canal in accordance with the clinical rules for an adhesive treatment, disinfect, dry and prepare for Total Etch-Technique (TE) or Self-Etch-Technique (SE).

Further on to Total Etch-Technique (3a. and 3b.) or Self-Etch-Technique (4.).

3 a. Total Etch-Technique

Apply FANTESTIC® ETCH (Order-No.: EGS3000) first on enamel (let it effect for about 20-30 s on enamel) and then on dentin (let it effect for about 10 s on dentin). In the case of teeth in high fluoride let it effect 90-120 s (or repeat the first step). Then rinse for about 20 s with water and dry the enamel with oil- and water-free air. The dentin should stay wet! Etched enamel should appear chalky white.

The root canal will be etched with FANTESTIC® ETCH and rinsed for about 20 s with water and then dried with paper tips.

3 b. Application of a dual-cure adhesive system like FANTESTIC® FLOWSIVE respectively - FLOWSIVE 2 each together with FANTESTIC® FLOWSIVE ACTIVATOR

Equivalent amounts of FANTESTIC® FLOWSIVE (Order-No.: FFF3100) respectively FANTESTIC® FLOWSIVE 2 (Order-No.: FFF3105) each together with FANTESTIC® FLOWSIVE ACTIVATOR (Order-No.: FAF3150) will be mixed for dual-cure light protected for about 15 s. Apply the mixture by using MICRO BRUSHES (Order-No.: MPB3981) onto the prepared dentin and enamel areas and agitate intensively for about 30 s and aerate for about 15 s in a weak air stream until no mobile liquid film is no longer visible. The mixture will auto-cure within 3 min (37 °C, 99 °F). It can be light-cured additionally for about 20 s. Repeat these steps for at least one time for achieving higher adhesive strength. Continue the adhesive treatment immediately.

4. Application of the self-etching, dual-cure adhesive system FANTESTIC® FLOWSIVE SE (Self-Etch)

FANTESTIC® FLOWSIVE SE (Order-No.: FFF5155) Liquid A and Liquid B will be mixed light protected 1:1 within the enclosed FANTESTIC® Mixing moulds (Order-No.: FFF5180) for about 10 s.

The homogenous mixture at this point is self etching (SE) for dentin and enamel.

Then the mixture will be applied with a brush onto the slightly wet dentin and onto the slightly wet enamel areas and will be agitated intensively for about 30 s (enamel) and about 15 s (dentin) and aerated for about 15 s in a weak air stream. The mixture can be light-cured (if possible) for about 20 s to get an optimal chemical adhesion. The self-cure system will auto-cure independently of light-curing within about 5 min (37 °C, 99 °F).

5 a. Preparation of FANTESTIC® Z CORE DC 5 ml-double-chamber syringe (KwikMix®)

Before first use screw off the original cap and waste it. Before screwing on a new mixing cannula the double-chamber syringe should be tested if both pastes will be extruded equally.

FANTESTIC® Z CORE DC (KwikMix®) will be mixed light protected and automatically using the MIXING CANNULA MCN 2 (Order-No.: MIK1041) + INTRA-ORAL-TIP MCN 2 (Order-No.: IOT2041) of the syringe system.

Instead of the INTRA-ORAL-TIP MCN 2 a fine root-canal application tip is available (ROOT-CANAL-TIP, Order-No.: RCT2041), with which the material could be applied directly and efficiently.

Further on instead of of MIXING CANNULA MCN 2 + INTRA-ORAL-TIP MCN 2 a NEEDLE MIXING CANNULA (Order-No.: MIK1041N) with an integrated, bendable application needle is available.

For this type of double-chamber syringe a mixing pistol for a comfortable application is available (MIXING PISTOL N 5, Order-No.: MIP1041).

5 b. Preparation of FANTESTIC® Z CORE DC 25 ml-cartridge (AutoMix)

Before first use screw off the original cap and waste it. Before screwing on a new mixing cannula the double chamber syringe should be tested if both pastes will be extruded equally.

FANTESTIC® Z CORE DC (AutoMix) will be mixed light protected and automatically using the yellow MIXING CANNULA LN (Order-No.: MIK1099) + INTRA-ORAL-TIP L (Order-No.: IOT2001) or using the blue MIXING CANNULA MLN (Order-No.: MIK2099) + INTRA-ORAL-TIP L (Order-No.: IOT2001).

Further on instead of of MIXING CANNULA LN + INTRA-ORAL-TIP L a yellow NEEDLE MIXING CANNULA LNN (Order-No.: MIK1099N) or a blue NEEDLE MIXING CANNULA MLNN (Order-No.: MIK2099N) with an integrated, bendable application needle is available.

For this 25 ml-cartridge a commercial mixing pistol (MIXING PISTOL KN, Order-No.: MIP1083) is required.

6. Pretreatments

6 a. Pretreatments of root canal posts

In front of an adhesive luting root canal posts must be silanized. Pay attention to the manufacturer's instructions.

6 b. Pretreatments of tooth structure (Total-Etch-Technique)

In this case a dual-cure adhesive system will be applied (like FANTESTIC® FLOWSIVE or FANTESTIC® FLOWSIVE 2 each with FANTESTIC® FLOWSIVE ACTIVATOR) onto the prepared dentin and enamel areas, see point 3 a. and 3 b).

6 c. Pretreatments of tooth structure (Self-Etch-Technique)

In this case a self-etching, dual-cure adhesive system like FANTESTIC® FLOWSIVE SE (Self-Etch) will be applied, see point 4.

Attention: If the FANTESTIC® Z CORE DC comes in contact with the adhesive system the chemical self-curing reaction will be accelerated.

7. Adhesive luting of root canal posts

For adhesive luting (cementing) of glass fiber posts or fiber reinforced composite (FANTESTIC® POST, Order-No.: FPS2X1X, x depends on geometry and size) the self-cure or dual-cure mode is valid.

The root canal will be prepared lege artis, disinfected (e. g. with 5 % NaOCl) and dried. Then the adhesive system according to section 3.-4. will be applied. Remove excess adhesive from the root canal using paper points.

Before application a small amount of material (size like rice corn to pea) should be pressed out and should be wasted immediately.

FANTESTIC® Z CORE DC will be applied onto the post (it can be silanized using FANTESTIC® CerBond, Order-No.: FCF4100 or ZM-Bond®, Order-No.: ZMF4200), and then the post will immediately be screwed into the root canal within the working time of about 1 min 30 s and hold by a slightly pressure of fingers.

FANTESTIC® Z CORE DC should not be luted directly in the root canal because of the applied adhesive system a premature polymerization could occur.

The self-cure system will auto-cure within about 3 min 30 s after working time. The coronal area can be light-cured optionally for about 40 s.

The filled MIXING CANNULA functions as a seal until the next application. Do not re-use the original cartridge cap.

Then the definite core build-up can be processed.

8. Definite core build-up

Before application a small amount of material (size like rice corn to pea) should be pressed out and should be wasted immediately.

FANTESTIC® Z CORE DC can be applied directly in thin layers within a working time of about 1 min 30 s (by using a matrix if necessary). The self-cure system will auto-cure within about 3 min 30 s after working time. Each layer placed can be light-cured optionally or in between for about 40 s.

In order to avoid the oxygen inhibition layer (about 100 µm) the water soluble FANTESTIC® OxyBlok (Order-No.: FOS3970) should be applied on the restoration margins immediately after removal of excess material. After final polymerization OxyBlok can be rinsed off with water.

The definite core preparation can be done using rotary instruments.

Precautions to application, product life and/or storage

Avoid contact of the etching gel with soft tissues, skin or eyes. If there is an accidental contact, flush immediately with numerous amounts of water.

Pay attention to the material safety data sheets (MSDS >> r-dental.com).

Do not use FANTESTIC® Z CORE DC in case of allergic reactions. Areas close to the pulp should be protected.

Avoid contact with soft tissues, skin or eyes. Do not adjust the viscosity by using any resin. Medical gloves do not protect against methacrylates.

Do not store FANTESTIC® Z CORE DC in proximity of eugenol containing materials nor let it come in contact with eugenol containing materials. Eugenol can impair hardening and cause discoloration. Further on no oxidic disinfection agents (like hydrogen peroxide) should be used, because these can inhibit hardening.

Attention: When using the dual-cure mode with FANTESTIC® Z CORE DC a self- or dual-cure adhesive system must be used. Do not use a self-etching, light-curing one-component adhesive system.

Before first use screw off the original cap and waste it. Before screwing on a new mixing cannula the double-chamber syringe should be tested if both pastes will be extruded equally.

Before application a small amount of material (size like rice corn to pea) should be pressed out and should be wasted immediately.

The filled MIXING CANNULA functions as a seal until the next application. Do not re-use the original cartridge cap.

Shelf life see label. Do not use after expiry date.

Before first use FANTESTIC® Z CORE DC must be stored within a refrigerator. After first use FANTESTIC® Z CORE DC should be stored lower 20 °C (68 °F) (optimal 8 °C - 20 °C, 46 °F - 68 °F) and horizontally. The opened FANTESTIC® Z CORE DC should be used up within three months, see lid.

Avoid exposure to sunlight. Keep out of the reach of children!

FANTESTIC, Cal de Luxe, ZM-Bond and KwikkMix are registered trademarks of R-dental Dentalerzeugnisse GmbH.

Caution: For dental use only. US Federal Law restricts this device to sale by or on the order of a dentist or licensed professional, to use or to order the use of this device.

The described product must be used according to the instructions for use and is determined for dental use only. There is no liability for defections resulting from inappropriate use or not intended purposes. The user is responsible for testing the suitability of the product for any purposes not explicitly stated in the instructions for use.

If there are damages to this product within the warranty period, R-dental is obliged to replace or to repair this product. There is no further claim. Outside of the liability by law there is no liability of R-dental for any loss or damages arising from this product, independently whether these damages are direct, indirect, special, incidental or consequential.

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© 2016 R-dental Dentalerzeugnisse GmbH. All rights reserved