R-SI-LINE® METAL-BITE®

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**Hersteller:**R-dental Dentalerzeugnisse GmbH • Winterhuder Weg 88 • 22085 Hamburg  
T 040-30707073-0 • F 040-30707073-73 • I www.r-dental.com**R-SI-LINE® METAL-BITE®**  
**REF MBK1074****Beschreibung und Indikation**

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R-SI-LINE® METAL-BITE® in automixing cartridges is an extremely fast, extremely hard and precise addition curing universal registration material on vinyl polysiloxane base (A-silicone), and is excellent applicable for bite-plates, too.

R-SI-LINE® METAL-BITE® is scanable and therefore usable for powderfree optical registration of data (antagonist registration), for example with CEREC®-cameras (SIRONA) (CAD/CAM/CIM-Impression Technology).

**Characteristics and composition**

R-SI-LINE® METAL-BITE® is thixotropic, very suitable for cutting and contouring, inflexible and not crumbly because of its high material lattice-constant. Registrations of R-SI-LINE® METAL-BITE® are permanently dimensionally stable.

R-SI-LINE® METAL-BITE® is delivered in auto-mixing cartridges and can be directly and bubblefree syringed onto the occlusal areas and so on.

**Recommended procedures****1. Preparing the auto-mixing cartridge/R-SI-LINE MIXING PISTOL N**

Introduce the auto-mixing cartridge into the MIXING PISTOL N (REF MIP1081) and lock it. Before first use screw off the original cap and dispose it. When screwing on a new mixing cannula the cartridge should be tested if both pastes will be extruded equally. Then screw on the mixing cannula R-SI-LINE MIXING CANNULA SXN (REF MIK1093). The mixing pistol is ready for application. For direct, broad application onto the teeth use the R-SI-LINE BITE-TIPS (REF BRT1000).

**2 a. Registration**

Before application a small amount of material (size of a pea) should be squeezed out and should be disposed immediately. The automatically mixed R-SI-LINE® METAL-BITE® is applied directly onto the occlusal areas or onto the bite-plate (picture 2 a.). The occlusion should be done within 20 s (icp). The registrate can be taken out of the mouth after 65 s.

**2 b. Data registration (static, icp-antagonist-registration)**

The automatically mixed R-SI-LINE® METAL-BITE® is applied directly onto the prepared occlusal areas. After hardening the icp-registrate will not be removed and will be prepared mesial and distal (picture 2 b.). Then the mesial, the registered (of the antagonist) and the distal occlusal area can be optically registered (Scan).

**2 c. Data registration (dynamic, FGP-antagonist-registration)**

After hardening the hikp-registrate will be prepared intraoral mesial and distal and will be reduced to a FGP-table with occlusal informations. Then fresh material will be applied onto this table and excursion movements (trusion movements) will be done. Then the dynamic registration can be scanned (literature >> www.r-dental.com).

Times mentioned for 74 °F. Higher temperatures shorten, lower temperatures prolong these times. The working time of the material can be prolonged by adjusting R-SI-LINE® METAL-BITE® down to the temperature of a refrigerator.

**3. Disinfection**

The registrates can be disinfected with standard disinfection solutions.

## Special storage and/or handling conditions

Before first use screw off the original cap and dispose it. Closing the cartridge by turning on the original cap between use could cause cross contamination of base - and catalyst paste and thus premature setting. When screwing on a new mixing cannula the cartridge should be tested if both pastes will be extruded equally.

Before application a small amount of material (size of a pea) should be squeezed out and should be disposed immediately.

The filled MIXING CANNULA functions as a seal until the next application. Do not re-use the original cartridge cap.

Do not touch tooth surfaces, preparations or retraction cords with latex gloves.

Expiry date see label. Do not use after expiry date.

Do not store above 78 °F!

Avoid exposure to sunlight. Keep out of the reach of children!

The registrates can be stored over months without shrinkage.

## Technical data

Mixing time:	auto-mixing
Working time incl. mixing time:	< 20 s
Recommended time in the mouth:	≥ 45 s
Time (start mixing - withdrawal of mouth):	≥ 65 s

R-SI-LINE and METAL-BITE are registered trademarks of R-dental Dentalerzeugnisse GmbH. CEREC is not a registered trademark of R-dental Dentalerzeugnisse GmbH.

**Caution:** For dental use only. US Federal Law restricts this device to sale by or on the order of a dentist or licensed professional, to use or to order the use of this device.

## Warranty

R-dental Dentalerzeugnisse GmbH warrants this product will be free from defects in material and manufacture. R-dental Dentalerzeugnisse GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. The user is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and R-dental Dentalerzeugnisse GmbH's sole obligation shall be repair or replacement of the product.

## Limitation of Liability

Except where prohibited by law, R-dental Dentalerzeugnisse GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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
## Material Safety Data Sheet according to Regulation (EU) 1907/2006/EC, Article 31 see www.r-dental.com

This short version: 05/2017

1. Identification of the substance/mixture and of the company/undertaking:  
Trade name: **R-SI-LINE ® METAL-BITE ®**  
**REF MBK1074**  
Relevant identified uses of the substance: Registration material  
Manufacturer/supplier: R-dental Dentalerzeugnisse GmbH - Winterhuder Weg 88 - 22085 Hamburg - DE.  
E-mail: info@r-dental.com  
Web: www.r-dental.com  
Emergency telephone number: Dr. A. Sprafke, T +49 40 - 30707073-0
2. Hazards identification:  
Classification of the substance: See chapter 15.
3. Composition/information of ingredients:  
Chemical characterization: Mixture of vinyl polysiloxanes/hydrogen-poly-dimethylsiloxanes/pigments/silica.  
Hazardous components: None.
4. First-aid measures:  
General information after skin contact: Wash with plenty of water and soap.  
General information after eye contact: Rinse with plenty of water and consult a physician.  
General information after swallowing: Seek medical advice immediately.
5. Fire fighting measures, Extinguishing media: CO<sub>2</sub>, water.
6. Accidental release measures:  
Personal precautions: Not necessary.  
Environmental precautions: Do not let material come into contact with the environment in a not cured state.  
Methods for cleaning up not cured material: Pick up the material.  
Methods for cleaning up cured material: Treat it like plastic disposal.
7. Handling and storage:  
Precautions for safe handling: None by correct handling, for dental use only, see instructions for use.  
Conditions for safe storage: Store horizontally at dry place, at room temperature, away from food stuffs.  
Requirements at store rooms and containers: Store dry at room temperature.  
Requirements for storage with other products: Avoid contact with acids, bases, alcohols.
8. Exposure controls/personal protection:  
The product contains no relevant quantities of components with limits of value to be supervised at the workplace:  
Personal protective equipment: Not necessary.  
General measures of protection and hygiene: Avoid contact with eyes, normal hygienic measures.  
Eyes: Protecting.  
Respiration: Not necessary.  
Hands: Carefully handling.
9. Physical and chemical properties:  
Appearance:  
Form: Paste.  
Colour: Grey.  
Odour: Odourless.  
Information on change in the physical state (Value, unit, method):  
Flash point: not applicable.  
Density: 1,60 g/cm<sup>3</sup> (20 °C, 68 °F).  
Viscosity: Heavy viscosity.  
pH: Neutral.  
Solubility in/miscibility with: A little bit soluble in toluene and petrolether.  
Water solubility: Not soluble.  
Content of solvents: None.  
Content of solids: Not applicable.
10. Stability and reactivity:  
Stability and reactivity: Stable.  
Incompatibility with other substances: Small evaluation of hydrogen traces with acids, bases, alcohols, powdered metals and metal oxides.  
Hazardous decomposition products: None by right use.
11. Toxicological information (estimated):  
Acute toxicity: Non toxic, no absorption.  
Primary irritation skin: Irritating (species rabbit).  
Eye: Irritating.  
Inhalation: Non volatile.  
Ingestion: Possible toxic when swallowed.  
Sensitization of product: Not observed.  
Additional toxicological information: During the long period of use as impression material no health diseases are reported.
12. Ecological information, general information:  
Degrance: Not determined.  
Fish toxicity: Not determined.  
Classification of water endangerment: Not determined.  
Additional information: Do not let the not cured product come into contact with the environment.
13. Disposal consideration:  
Recommendation: Not cured material use proper landfill disposal in accordance with local, state and federal regulations.
14. Transport information:  
Transport information: Not classified as dangerous goods. Keep away from stuff, acids and bases.  
Land transport RID/ADR, GGVS/GGVE: Not necessary.  
Ship transport ADN/ADNR: Not necessary.  
Sea transport IMDG/GGVSea: Not necessary.  
Air transport ICAO-TI and IATA-DGR: Not necessary.
15. Regulatory information, classification according to EEC-guidelines:  
The product is a class I Medical Device according to the MDD 93/42/EEC.  
According to GEFStoffV/EG: No classification necessary to actual EEC guidelines.
16. Other informations:  
The above information is based on our present day knowledge and relates solely to the safety requirements of the product. The data do not signify any warranty with regards to product properties. However users of the product should satisfy themselves that the information given above is sufficient and correct for their circumstances of use.

Prepared by: Dr. Andreas Sprafke, R-dental Dentalerzeugnisse GmbH.

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© 2021 R-dental Dentalerzeugnisse GmbH. All rights reserved.