FANTESTIC® FLOWFILL SE

CE 0482

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R-dental Dentalmaterialien GmbH · Winterhuder Weg 88 · 22085 Hamburg  
T 040-30707073-0 · E info@r-dental.com · I www.r-dental.com

## FANTESTIC® FLOWFILL SE

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### FANTESTIC® FLOWFILL SE

REF FFS31XX (in syringes, XX depends on shade)

#### Product description

FANTESTIC® FLOWFILL SE is a flowable, self-adhesive and light-cure composite for anterior and limited posterior restorations and for sealings. Beside its good mechanical properties, the physical properties of FANTESTIC® FLOWFILL SE are comparable to conventional flowable composites. The optimized flow behavior enables an easy application. The self-adhesive properties of FANTESTIC® FLOWFILL SE require no etching, priming and bonding of enamel and dentin in the restorative therapy. Only for an application as a pit and fissure sealant, enamel etching for cleaning is recommended.

FANTESTIC® FLOWFILL SE is based on methacrylates and inorganic filler particles of a size of 0.05-1.0 µm. The total filler load is 59 wt.-% and 45 vol.-%. The radiopacity of FANTESTIC® FLOWFILL SE is 150 % aluminium (radiopacity of enamel: ~ 200 % aluminium, of dentin: ~ 100 % aluminium).

FANTESTIC® FLOWFILL SE meets the requirements of DIN EN ISO 4049, type 1, class 2, group 1.

#### Indications/intended use

FANTESTIC® FLOWFILL SE is indicated for

- small restorations of class I,
- lining of class I and class II restorations,
- pit & fissure sealings,
- extended pit & fissure sealings and
- filling of undercuts.

#### Performance features

The performance features of the product meet the requirements of the intended use.

#### Contraindications

Irritations resulting from direct contact with the pulp are possible. Therefore for pulp protection areas close to the pulp must be covered with a thin layer of pulp capping material, see section 4.

#### Patient target group

Persons who are treated during a dental procedure.

#### Intended users

This medical device should only be used by a professionally trained dental practitioner.

#### Incompatibilities with other materials

Do not use FANTESTIC® FLOWFILL SE in combination with eugenol containing materials. Eugenol impairs the polymerization of the composite. Do not store FANTESTIC® FLOWFILL SE in proximity of eugenol containing materials. Unpolymerized FANTESTIC® FLOWFILL SE should not come in contact with eugenol containing materials.

#### Application

##### 1. Preparation

Prior to the preparation a cleaning of the tooth with a fluoride-free paste (e. g. fluor of pumice with water) is recommended.

##### 2. Shade matching

Ascertain the tooth shade while teeth are still moist and select the appropriate FANTESTIC® FLOWFILL SE shade (FANTESTIC® FLOWFILL SE-products are available in VITA®-shades).

##### 3. Isolation

A dry area and use of rubber dam is recommended.

##### 4. Preparation of the cavity

The cavity will then be prepared in accordance with clinical rules and with minimal tooth reduction. If possible, margins should have a slight (0.5-1.0 mm) bevel placed in the enamel to increase the surface area for higher bond strength.

For pulp protection areas close to the pulp should be protected with a thin layer of calcium hydroxide, like Cal de Luxe®, according to the instructions for use or with a bioactive pulp capping material, like ReViCal® (REF RVC100X) or ReViCal® Liner (REF RVC10XX), according to the instructions for use.

##### 5. Application of the flowable composite FANTESTIC® FLOWFILL SE

Screw off the original cap of the syringe and substitute it by one of the supplied application cannula (REF FFK3281) by turning clockwise. Store the screw original-cap. Squeeze out the first small drop on a pad to ensure safe handling.

After use the application cannula will be removed and again the original-cap replaced.

The application cannulas are for single use only.

##### 5 a. Small class I restorations and lining of class I and II restorations

FANTESTIC® FLOWFILL SE will be applied with the application cannula in a thin layer onto the prepared areas. By using a MICRO BRUSH (REF MPB3981) FANTESTIC® FLOWFILL SE will be worked onto the entire cavity walls and beveled areas under slight pressure for 20 s. Excess material will be removed around the margins with a microbrush.

FANTESTIC® FLOWFILL SE will be light-cured for **20 s** for lighter shades (e. g. A2) and **30 s** for darker shades (e. g. A3).

After lining of the cavity walls and the beveled areas the restoration will be build up with FANTESTIC® FLOWFILL SE in layers of max. 2 mm.

Light-cure each increment separately. FANTESTIC® FLOWFILL SE will be light-cured for **20 s** for lighter shades (e. g. A2) and **30 s** for darker shades (e. g. A3).

##### 5 b. Pit and fissure sealing

Before sealing the enamel will be cleaned with a fluoride-free paste and rinsed thoroughly, dried in a water and oil-free air stream for 5 s leaving a slightly moist shining surface on the enamel.

Note: In addition to the cleaning as described above, etching of enamel with a 37 % phosphoric acid etching gel is recommended. Etch the enamel surface for 15-20 s, rinse for 10 s and dry in a water and oil-free air stream.

FANTESTIC® FLOWFILL SE will be placed with an application cannula in a thin layer directly onto the prepared areas. By using a MICRO BRUSH (REF MPB3981) FANTESTIC® FLOWFILL SE will be worked onto the prepared areas under slight pressure for 20 s. Excess material will be removed around the margins with a micro brush.

FANTESTIC® FLOWFILL SE can be light-cured for **20 s** for lighter shades (e. g. A2) and **30 s** for darker shades (e. g. A3). Hold the light emission window as close as possible to the filling material. The occlusion will be checked and corrected, if necessary.

##### 6. Finishing and high gloss polishing

Depending on the shape of the area contour FANTESTIC® FLOWFILL SE can be contoured with finishing diamonds, flexible grinding discs or carbide burs. Polish to high gloss with polishing discs or polishing brushes.

##### 8. Storage

Do not store above 25 °C (77 °F)! Protect from direct sunlight. Close syringe immediately after use to avoid exposure to light.

Do not use after expiry date.

##### Additional Notes/Warnings

FANTESTIC® FLOWFILL SE can be light-cured with commercial polymerization units (halogen or LED) with a wave length of 400-500 nm and with a light intensity of at least 1,000 mW/cm².

The ambient light of the dental lamp may start polymerization of the composite.

Do not adjust the viscosity of the composite restorative material by using any resin.

Contact of resin pastes with skin should be avoided, especially by anyone having known resin allergies.

Contact with eyes may cause severe eye damage. Wear eye protection. In case of contact with eyes rinse immediately with plenty of water and seek medical advice.

Commercial medical gloves do not protect against sensitizing effect of methacrylates.

\*VITA is a registered trademark of VITA-Zahnfabrik, Bad Säckingen.

Keep out of the reach of children!

##### Composition

Dimethacrylates, dental glass, silicone dioxide, photo initiators.

##### Disposal

Disposal of the product according to local authority regulations.

##### Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

##### Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

FANTESTIC, Cal de Luxe and ReViCal are registered trademarks of R-dental Dentalerzeugnisse GmbH.

##### Warranty

R-dental Dentalmaterialien GmbH warrants this product will be free from defects in material and manufacture. R-dental Dentalmaterialien GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. The user is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and R-dental Dentalmaterialien GmbH's sole obligation shall be repair or replacement of the product.

##### Limitation of Liability

Except where prohibited by law, R-dental Dentalmaterialien GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

© 2024 R-dental Dentalmaterialien GmbH. All rights reserved.