FANTESTIC® FLOWFILLCEM DC

CE 0482

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Hersteller:R-dental Dentalmaterialien GmbH · Winterhuder Weg 88 · 22085 Hamburg
T 040-30707073-0 · F 040-30707073-73 · I www.r-dental.com**FANTESTIC® FLOWFILLCEM DC (KwikkMix®)**

REF FFS3266 (A3 Universal)

REF FFS3267 (weiß-opak)

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REF FFS3267 (white-opaque)

Description and indication

FANTESTIC® FLOWFILLCEM DC is a flowable, self-curing polymer based dental luting material (micro-hybrid composite) with light-curing option (dual-cure, DC) for definite intraoral adhesive luting ("cementation") of indirect restorations like inlays, onlays, crowns, bridges and veneers.

Characteristics and composition

FANTESTIC® FLOWFILLCEM DC is very radiopaque with a very high strength and a relatively high filler content. It is highly polishable, and it is very esthetic with an excellent colour stability. It is based on methacrylates and inorganic filler particles of 0.05-1.0 µm. The filler content is about 61 % by weight and 41 % by volume. The radiopacity of FANTESTIC® FLOWFILLCEM DC is 150 % aluminium (radiopacity of enamel: ~ 200 % aluminium, of dentin: ~ 100 % aluminium).

FANTESTIC® FLOWFILL meets the requirements of DIN EN ISO 4049, type 2, class 2, group 1.

Recommended procedures

1. Isolation

A dry area and use of rubber dam (REF KDL31XX, XX = depends on strength) or latex-free dam (REF KDL3500) is recommended.

2. Preparation of 5 ml-double-chamber syringe (KwikMix®)

Before first use screw off the original cap and waste it. Before screwing on a new mixing cannula the double-chamber syringe should be tested if both pastes will be extruded equally. Then the mixing cannula is put on and locked by turning it 90 ° clockwise.

FANTESTIC® FLOWFILLCEM DC (KwikMix®) will be mixed light protected and automatically using the MIXING CANNULA MCN 2 (REF MIK1041) + INTRA-ORAL-TIP MCN 2 (REF IOT2041) of the syringe system.

Instead of the INTRA-ORAL-TIP MCN 2 a fine root-canal application tip is available (ROOT-CANAL-TIP, REF RCT2041), with which the material could be applied directly and efficiently.

Further on instead of of MIXING CANNULA MCN 2 + INTRA-ORAL-TIP MCN 2 a NEEDLE MIXING CANNULA (REF MIK1041N) with an integrated, bendable application needle is available.

For this type of double-chamber syringe a mixing pistol for a comfortable application is available (REF MIP1041).

3. Pretreatments

3 a. Pretreatment of luting areas of indirect restorations

In front of an adhesive luting glass-ceramic restorations must be silanized according to the instructions for use (FANTESTIC® CerBond, REF FCF4100).

Restorations made of zirconium dioxide, metal (precious, non-precious) and metal-supported restorations must be treated in front of an adhesive luting with the universal primer ZM-Bond® (REF ZMF4200). ZM-Bond® is a single-component primer, which promotes a chemical adhesion between luting composites and restorative materials.

Attention: Do not use phosphoric acid on the adhesive contact areas of zirconium dioxide, metal (precious, non-precious), metal-supported restorations!

3 b. Pretreatments of tooth structure

On teeth FANTESTIC® FLOWFILLCEM DC must applied together with dual-cure or self-etching, dual-cure adhesive systems according to the instructions for use (for example FANTESTIC® FLOWSIVE (REF FFF3100) or - FLOWSIVE 2 (REF FFF3105) with - ACTIVATOR (REF FAF3150) or FANTESTIC® FLOWSIVE SE (REF FFF5155)).

4. Adhesive luting of inlays, onlays, crowns and bridges

FANTESTIC® FLOWFILLCEM DC will then be mixed automatically and light protected using the mixing cannulas of the safety cartridge system. Before every application a small amount of material (2-3 mm, size like pepper corn) should be pressed out and should be wasted immediately.

FANTESTIC® FLOWFILLCEM DC will be applied directly in a thin layer of about 0.5 mm onto the primed areas of the restoration and of the primed teeth and place the restoration with slight pressure within 2 min (working time incl. mixing time at 74 °F (23 °C) and self-curing) and hold until it is set completely. Excess material of FANTESTIC® FLOWFILLCEM DC can be removed by using brushes, FANTESTIC® MICRO BRUSH (REF MPB3981), dental floss or scaler.

After removal of excess material the marginal areas will be light-cured with a polymerization lamp (wave length range 400-500 nm) with a light intensity of at least 1,000 mW/cm² from each direction for 20 s.

FANTESTIC® FLOWFILLCEM DC will polymerize completely within about 4 min, independently of light-curing.

In order to avoid the oxygen inhibition layer (about 100 µm) the water soluble FANTESTIC® OxyBlock (REF FOS3970) should be applied according to the instructions for use on the restoration margins immediately after removal of excess material. After final polymerization OxyBlock can be rinsed off with water.

Then the margins should be polished.

4. Adhesive luting of veneers

The veneers will be adapted with try-in gel to the prepared areas of the tooth. Avoid gingival contact to prevent bleeding. After removal of the veneers the prepared enamel areas and the veneers will be rinsed intensively with water and dried with oil-free air.

The tooth which should be veneered will be isolated with interproximal strips in order to protect adjacent teeth (not being veneered) from etching material and adhesive.

Apply an etching gel (e. g. FANTESTIC® ETCH, REF EGS3000) onto the enamel areas and let it effect for 15 s and rinse intensively with water. After drying place a thin coat of a light-cure adhesive (e. g. FANTESTIC® Universal Bond, REF FUB6000) according to the instructions for use onto the etched enamel areas. Then FANTESTIC® FLOWFILLCEM DC will be applied directly.

The interproximal strips will be removed prior to the application of the veneers. This step assures a complete and automatic setting of the veneers.

After placing the veneer light-cure with a polymerization lamp (wave length range 400-500 nm) with a light intensity of at least 1,000 mW/cm² through the tooth structure and the veneer for 20 s.

In cases the tooth structure and the veneer are too opaque and the light-curing is not sufficient for polymerization, FANTESTIC® FLOWFILLCEM DC self-cures within 4 min.

After finishing polymerization and removal of excess material the margins should be polished.

Precautions to application, product life and/or storage

Avoid contact of the etching gel with soft tissues, skin or eyes. If there is an accidental contact, flush immediately with numerous amounts of water.

Do not use FANTESTIC® FLOWFILLCEM DC in case of allergic reactions to methacrylates. Avoid contact of FANTESTIC® FLOWFILLCEM DC with soft tissues, skin or eyes. Do not adjust the viscosity by using any resin.

Do not store FANTESTIC® FLOWFILLCEM DC in proximity of eugenol containing materials nor let it come in contact with eugenol containing materials. Eugenol can impair hardening and cause discoloration.

Before first use screw off the original cap from the double-chamber syringe and waste it. Before screwing on a new mixing cannula the double-chamber syringe should be tested if both pastes will be extruded equally.

Before every application a small amount of material (2-3 mm, size like pepper corn) should be pressed out and should be wasted immediately.

The filled MIXING CANNULA functions as a seal until the next application. Do not re-use the cartridge cap.

Pay attention to the material safety data sheets (MSDS >> r-dental.com).

Shelf life see label. Do not use after expiry date.

Do not store the double-chamber syringe above 68 °F (20 °C) and store it horizontally, that means lying. Store unopened material in the refrigerator. Opened FANTESTIC® FLOWFILLCEM DC should be used within 3 months.

Avoid exposure to sunlight. Keep out of the reach of children!

FANTESTIC, ZM-Bond and KwikMix are registered trademarks of R-dental Dentalerzeugnisse GmbH.

Caution: For dental use only.

Warranty

R-dental Dentalmaterialien GmbH warrants this product will be free from defects in material and manufacture. R-dental Dentalmaterialien GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. The user is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and R-dental Dentalmaterialien GmbH's sole obligation shall be repair or replacement of the product.

Limitation of Liability

Except where prohibited by law, R-dental Dentalmaterialien GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

© 2021 R-dental Dentalmaterialien GmbH. All rights reserved.