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FANTESTIC® POST

CE 0546

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Instructions for use

English - FPS2X1X - Rev. 7 - 2021/09

1) After root canal filling, take care not to use cement containing eugenol, you can prepare the canal with traditional instruments. We suggest to use a rubber dam.

2) Once you have chosen the post diameter, you can proceed to bore using the calibrated drill for a correct post positioning. Drive the post 6 - 7 mm deep in the root canal. The post should be placed without any effort and oscillations.

Specific recommendation for posts of Ø 0,8 mm.: these models can be used, on discretion of the dentists, only in conditions of limited stress (e.g. in incisors or in connection with other posts).

3) Once you have determined the correct length of the post, cut it with a diamond disk. Never use drills or nippers. It is recommended to carry out the cut out of the oral cavity. It is suggested to wash the device after the cut in order to eliminate possible dust residues. We recommend to clean the device with alcohol before use or in case of accidental or improper contact.

4) After the mechanical treatment proceed with the chemical treatment. We suggest to etch the dentin with the combined action of EDTA and sodium hypochlorite for a correct root canal irrigation. Rinse and dry softly with paper cones and / or air jets. Do not dehydrate.

5) Handle the post with suitable instruments only in order to avoid the hands contact. In this case, clean the post with alcohol or with a degreasing agent. Coat the canal walls with dentin-enamel adhesive. Use a cement-composite with a resin base for affinity with the resin of the post. The consistency of the paste should be fluid in order to properly bind to the dentinal tubules. Insert the cement into the canal using a Lentulo. Light-cure the visual part for 20 seconds and wait until self-curing is completed (about 4 minutes).

6) Core build-up can be performed on the post using light and self-curing composites.

7) The post cannot come into direct contact with the patient's tissues and does not replace the implant abutment.

8) Composition: Resin epoxy matrix, glass fibers and barium sulphate, carbon fibers, glass translucent fibers.

The device must be used only by professional dentists.

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