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Deutsch - RVC10XX-01/23

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ReViCal® Liner

CE 0482

Hersteller:

R-dental Dentalmaterialien GmbH • Winterhuder Weg 88 • 22085 Hamburg
T 040-30707073-0 • F 040-30707073-73 • I www.r-dental.com

ReViCal® Liner

REF RVC1005 (Einzelpackung, 1 Spritze RVC1005)

REF RVC1010 (Doppelpackung, 2 Spritzen RVC1005)

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REF RVC1010 (Double pack, 2 Syringes RVC1005)

Description and Characteristics

ReViCal® Liner is a bioactive, alkaline, light-cure, reinforced liner with MTA-fillers.

ReViCal® Liner can be applied very precisely with the delivered application cannulas. After that a filling can be placed immediately.

Indications/Intended use

ReViCal® Liner is indicated

1. as a liner (thin layer underfilling) and
2. for indirect pulp capping (Cp).

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

In singular cases, the material may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

Patient target group

Persons who are treated during a dental procedure.

Intended users

This medical device should only be used by a professionally trained dental practitioner.

Recommended procedures

1. Isolation

A dry area and use of rubber dam (REF KDL31XX, XX = depends on strength) or latex-free dam (REF KDL3500) is recommended.

2. Lege artis cavity preparation

The cavity will be prepared lege artis.

3. Application of ReViCal® Liner

When used as a liner (thin underfilling) and for indirect pulp capping apply first a light-cure dentin-enamel-adhesive according to the instructions for use of the manufacturer to ensure a fast and an effective adhesion.

First screw off the cap of the syringe and substitute it by an application cannula by turning 90° clockwise.

Store the screw cap.

The first small drop will be squeezed out on a pad to ensure a safe handling of the ReViCal® Liner.

Then ReViCal® Liner will be applied in one layer of ≤ 1 mm.

ReViCal® Liner should not be applied to the cavity edges and remove any paste sticking to them.

Due to the high opacity ReViCal® Liner will be light-cured for 40 s with a commercial polymerization unit (wave length range 400-500 nm) with a light intensity of at least 1,000 mW/cm².

After use the application cannula of the ReViCal® Liner-syringe will be removed and again the screw cap replaced.

For hygienic reasons the application cannulas are for single use only.

Storage

Store at 2 °C - 8 °C (36 °F - 46 °F) within the original seal-pack. After opening the original seal-pack waste it. Close the ReViCal® Liner-syringe immediately after first use and store it protected from light and moisture at room temperature 19 °C - 25 °C (66 °F - 77 °F) without original seal-pack. Once opened, ReViCal® Liner must be used within 6 months and before the expiry date. Shelf life see label. Do not use after expiry date. Keep out of the reach of children!

Precautions to application/warnings

ReViCal® Liner contains polymerizable methacrylate monomers. Unpolymerized ReViCal® Liner may have an irritating effect and may lead to a sensitizing reaction against methacrylates.

In singular cases, the material may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

Avoid contact of ReViCal® Liner with mucous membrane, skin or eyes. In the case of skin contact wash immediately with water. In case of contact with eyes, rinse immediately with copious amounts of water and consult a physician.

Do not use ReViCal® Liner in combination with substances containing eugenol. Eugenol inhibits the polymerization of the material. Neither store ReViCal® Liner in proximity of open eugenol containing products. Do not let unpolymerized ReViCal® Liner allow coming into contact with materials containing eugenol.

Commercial medical gloves do not protect against a sensitizing effect of methacrylates.

Pay attention to the material safety data sheets (MSDS >> r-dental.com).

ReViCal is a registered trademark of R-dental Dentalerzeugnisse GmbH.

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Warranty

R-dental Dentalmaterialien GmbH warrants this product will be free from defects in material and manufacture. R-dental Dentalmaterialien GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. The user is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and R-dental Dentalmaterialien GmbH's sole obligation shall be repair or replacement of the product.

Limitation of Liability

Except where prohibited by law, R-dental Dentalmaterialien GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

© 2023 R-dental Dentalmaterialien GmbH. All rights reserved.