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ArtOral® Chroma BLOKK

CE 0482

**Hersteller:**R-dental Dentalerzeugnisse GmbH • Winterhuder Weg 88 • 22085 Hamburg
T 040-30707073-0 • F 040-30707073-73 • I www.r-dental.com**ArtOral® Chroma BLOKK**

REF ACS1005 (Spritze)

REF ACT1005 (Tip) → ⊗

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REF ACT1005 (tip)

Description and indication

ArtOral® Chroma BLOKK is a light-cure, radioopaque restoration material (composite) that is used in combination with ArtOral® Chroma FILL to mask discolored teeth.

Characteristics and composition

ArtOral® Chroma BLOKK contains dimethacrylate based resins and inorganic filler with a particle size < 0.2 µm. The total filler content is 80 wt.-% and 76 vol.-%.

The radiopacity of ArtOral® Chroma BLOKK is 200 % Aluminium (radiopacity enamel ≈ 200 % Aluminium, dentin ≈ 100 % Aluminium).

The material properties of ArtOral® Chroma BLOKK meet the requirements of ISO 4049, type 1, class 2, group 1.

Indications

ArtOral® Chroma BLOKK is indicated

- for masking an oral cavity in case of extended class III and IV cavities and
- for masking of tooth discoloration.

Side effects

Irritations resulting from direct contact with the pulp are possible. Therefore for pulp protection areas close to the pulp must be covered with a thin layer of calcium hydroxide or a bioactive pulp capping material, see section 3.

Incompatibilities

Do not use ArtOral® Chroma BLOKK and ArtOral® Chroma FILL in combination with eugenol containing materials. Eugenol impairs the polymerization of the composite. Do not store ArtOral® Chroma BLOKK and ArtOral® Chroma FILL in proximity of eugenol containing materials. Unpolymerized ArtOral® Chroma BLOKK and ArtOral® Chroma FILL should not come in contact with eugenol containing materials.

Recommended procedures

1. Preparation of the tooth

It is recommended to clean the tooth with a fluoride-free polishing paste (e. g. fluor of pumice with water).

2. Isolation

A dry area and use of rubber dam (REF KDL31XX, XX = depends on strength) or latex-free dam (REF KDL3500) is recommended.

3. Preparation of the cavity

The cavity will then be prepared in accordance with clinical rules and with minimal tooth reduction.

For pulp protection areas close to the pulp should be protected with a thin layer of calcium hydroxide, Cal de Luxe®, (Ca(OH)₂, "ready-to-use" in a portion blister, REF CDL2503) according to the instructions for use or with a pulp cap material, a bioactive pulp capping material, like ReViCal (REF RVC1002) according to the instructions for use.

Then apply the Total-Etch-Technique (Etch & Rinse) or the Self-Etch-Technique with common adhesives for dentin and enamel according to the instructions for use.

The application of FANTESTIC® UniversalBond (REF FUB6000) as a self-etching one-component adhesive for use for Self-Etch-Technique or the selective enamel etching in combination with an etching gel according to the instructions for use is recommended.

4. Application of the composite ArtOral® Chroma BLOKK

In case of extended class III and class IV palatal cavities use ArtOral® Chroma BLOKK to ensure that the restoration may not appear dark. ArtOral® Chroma BLOKK can also be used as a first layer for class I, II and V cavities of dark or discolored teeth.

In case of repair of ceramic or composite restorations use a light-cure opaquer as first layer in case of chipped metal or zirconia surfaces.

Apply ArtOral® Chroma BLOKK in the cavity with a suitable instrument in layers of max. 1 mm. Light-cure each increment separately. ArtOral® Chroma BLOKK is easy mouldable. The use of a matrix is recommended.

Light-cure according to chapter 6.

5. Application of the composite ArtOral® Chroma FILL

The restoration will then be finished with ArtOral® Chroma FILL (REF ACS1000 syringe or REF ACT1000 tip according to the instructions for use.

6. Light-cure

ArtOral® Chroma BLOKK und ArtOral® Chroma FILL can be light-cured with commercial polymerization units (halogen or LED) with a wave length of 400-500 nm and with a light intensity > 1000mW/cm². Cure each layer for **20 s**. Hold the light emission window as close as possible to the filling material.

7. Storage

Do not store above 77 °F! Avoid exposure to sunlight.

Precautions to application, product life and/or storage

Avoid contact of calcium hydroxide Cal de Luxe® and the pulp capping material ReViCal with soft tissues, skin or eyes. If there is an accidental contact, flush immediately with copious amounts of water.

Pay attention to the Material Safety Data Sheets (MSDS >> www.r-dental.com).

The ambient light of the dental lamp may start polymerization of the composite.

Do not adjust the viscosity of the composite restorative material by using any resin.

Avoid contact with skin, mucous membrane and eyes.

Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates.

Commercial medical gloves do not protect against methacrylates.

Close syringe immediately after use to avoid exposure to light. Composites in tips are for single use only.

Expiry date see label. Do not use after expiry date.

ArtOral is a registered trademark of R-dental Dentalerzeugnisse GmbH.

Keep out of the reach of children!

Caution: For dental use only.

Warranty

R-dental Dentalerzeugnisse GmbH warrants this product will be free from defects in material and manufacture. R-dental Dentalerzeugnisse GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. The user is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and R-dental Dentalerzeugnisse GmbH's sole obligation shall be repair or replacement of the product.

Limitation of Liability

Except where prohibited by law, R-dental Dentalerzeugnisse GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

© 2021 R-dental Dentalerzeugnisse GmbH. All rights reserved.