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ReViCal®

CE 0482

Hersteller:

R-dental Dentalmaterialien GmbH · Winterhuder Weg 88 · 22085 Hamburg
T 040-30707073-0 · F 040-30707073-73 · I www.r-dental.com

ReViCal®

REF RVC1001 (Einzelpackung, 1 Spritze RVC1001)

REF RVC1002 (Doppelpackung, 2 Spritzen RVC1001)

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REF RVC1002 (double pack, 2 Syringes RVC1001)

Description and Characteristics

ReViCal® is a bioactive, light-cure, reinforced pulp capping material with MTA-fillers. It is designed to perform as a barrier and to protect the pulp. ReViCal® is thixotropic and can be applied very precisely with the delivered application cannulas even in deep cavity preparations. Because of light-cure a filling can be light-cured under control and can be placed immediately.

Indications/Intended use

ReViCal® is indicated for

1. direct pulp capping for any pulpal exposures, including after a

- carious pulp exposure,
- mechanical pulp exposure and
- a pulp exposure due to a trauma.

2. indirect pulp capping agent in deep preparations under

- amalgam restorations,
- class I and II composite restorations,
- cements and
- as an alternative to calcium hydroxide.

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

In singular cases, the material may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

Patient target group

Persons who are treated during a dental procedure.

Intended users

This medical device should only be used by a professionally trained dental practitioner.

Precautions

Because of its high pH value ReViCal® should not be used in direct intraoral environment. ReViCal® should not be placed on enamel or cavity margins or leave exposed material to intraoral environment.

The success of pulp capping procedures may be affected by preparation contamination. Use of a rubber dam or adequate isolation is recommended.

Pulpal hemorrhage must be controlled prior to the application of ReViCal®. If the hemorrhage cannot be controlled, an endodontic therapy must be considered.

Recommended procedures

First screw off the cap of the syringe and substitute it by an application cannula. Store the screw cap.

A. Direct pulp capping for any pulpal exposures

1. Isolation

A dry area and use of rubber dam (REF KDL31XX, XX = depends on strength) or latex-free dam (REF KDL3500) is recommended.

2. Lege artis cavity preparation

The cavity will be prepared lege artis.

When opening the pulp, first stop any bleeding that may occur. Then remove excess moisture with a sterile cotton pellet. The surface must be visibly moist.

3. Application of the pulp capping material ReViCal®

First screw off the cap of the syringe and substitute it by an application cannula. Store the screw cap.

Then ReViCal® will be applied punctiform and small-area with a layer of ≤ 1 mm to the opened pulp. All exposed areas of the opened pulp will be covered with ReViCal® and extended optionally up to 1 mm into the dentin.

Due to the high opacity ReViCal® will be light-cured for 40 s with a commercial polymerization unit (wave length range 400-500 nm) with a light intensity of at least 1,000 mW/cm².

Then a suitable, light-cure adhesive will be applied onto the cured ReViCal® layer according to the instructions for use of the manufacturer. Then it can be proceed with the restoration.

Note: When using an Etch & Rinse adhesive take care to etch only the exposed cavity areas. Then apply the adhesive to all cavity areas.

B. Indirect pulp capping in deep preparations

1. Isolation

A dry area and use of rubber dam (REF KDL31XX, XX = depends on strength) or latex-free dam (REF KDL3500) is recommended.

2. Lege artis cavity preparation

The cavity will be prepared lege artis.

Then all infected carious tooth structure must be removed. The surfaces must be visibly moist.

3. Application of the pulp capping material ReViCal®

First screw off the cap of the syringe and substitute it by an application cannula. Store the screw cap.

Then ReViCal® will be applied in small-area with a layer thickness of ≤ 1 mm directly onto the area close to the pulp.

Due to the high opacity ReViCal® will be light-cured with a commercial polymerization unit (wave length range 400-500 nm) with a light intensity of at least 1,000 mW/cm² for 40 s.

Then a suitable, light-cure adhesive will be applied onto the cured ReViCal® layer according to the instructions for use of the manufacturer. Then it can be proceed with the restoration.

Note: When use an Etch & Rinse adhesive take care to etch only the exposed cavity areas. Then apply the adhesive to all cavity areas.

After use the application cannula of the ReViCal® syringe will be removed and again the screw cap replaced.

The application cannulas are for single use only.

Storage

Store at 2 °C - 8 °C (36 °F - 46 °F) within the original seal-pack. After opening the original seal-pack waste it. Close the ReViCal®-syringe immediately after first use and store it protected from light and moisture at room temperature 19 °C - 25 °C (66 °F - 77 °F) without original seal-pack. Once opened, ReViCal® must be used within 6 months and before the expiry date. Shelf life see label. Do not use after expiry date. Keep out of the reach of children!

Precautions to application/warnings

ReViCal® contains polymerizable methacrylate monomers. Unpolymerized ReViCal® may have an irritating effect and may lead to a sensitizing reaction against methacrylates.

Do not use with patients with allergic reactions against acrylates.

Avoid contact of ReViCal® with mucous membrane, skin or eyes. In the case of skin contact wash immediately with water. In case of contact with eyes, rinse immediately with copious amounts of water and consult a physician.

Do not use ReViCal® in combination with substances containing eugenol. Eugenol inhibits the polymerization of the material. Neither store ReViCal® in proximity of open eugenol containing products. Do not let unpolymerized ReViCal® allow coming into contact with materials containing eugenol.

Commercial medical gloves do not protect against a sensitizing effect of methacrylates.

Pay attention to the material safety data sheets (MSDS >> r-dental.com).

ReViCal is a registered trademark of R-dental Dentalerzeugnisse GmbH.

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Warranty

R-dental Dentalmaterialien GmbH warrants this product will be free from defects in material and manufacture. R-dental Dentalmaterialien GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. The user is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and R-dental Dentalmaterialien GmbH's sole obligation shall be repair or replacement of the product.

Limitation of Liability

Except where prohibited by law, R-dental Dentalmaterialien GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

© 2023 R-dental Dentalmaterialien GmbH. All rights reserved.