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R-SI-LINE® MEDIUM SH

CE 0482

R-dental Dentalmaterialien GmbH · Winterhuder Weg 88 · 22085 Hamburg  
T 040-30707073-0 · E info@r-dental.com · I www.r-dental.com

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REF MEK1030SH (R-SI-LINE MEDIUM SH (super hydro))

REF MEK1030SHFS (R-SI-LINE MEDIUM SH FS (super hydro, fast set))

### Product description

R-SI-LINE® MEDIUM SH (- FS) in auto-mixing cartridges is used particularly as a correction material for the sandwich-technique (1-step/2-phase) in conjunction with the base impression material R-SI-LINE® PUTTYSOFT (REF PUS1020), R-SI-LINE® PUTTYSOFT FS (REF PUS1020FS), R-SI-LINE® PUTTY MATIC (REF PMK1013), R-SI-LINE® HEAVY MATIC (REF HMK1053) or R-SI-LINE® HEAVY KINGSIZE (REF HKK1051).

R-SI-LINE® MEDIUM SH (- FS) is a light to medium viscosity addition curing precision impression material on vinylpolysiloxane base. It can be directly and bubble-free applied onto the teeth or the tray.

It is very thixotropic and shows a controlled flow behavior. Thus it stays in the position on the tray, but under pressure it flows.

R-SI-LINE® MEDIUM SH (super hydro) and R-SI-LINE® MEDIUM SH FS (super hydro, fast set) are very hydrophilic and show a finest detail reproduction and an excellent detail reproduction of subgingival areas and other moist surfaces.

R-SI-LINE® MEDIUM SH (- FS) meet the requirements ISO 4823, Type 3 (light-bodied consistency).

### Indications/intended use

R-SI-LINE® MEDIUM SH (- FS) is indicated for

- one-step (1-step/2-phase) sandwich-technique with R-SI-LINE® PUTTYSOFT or R-SI-LINE® PUTTYSOFT FS (s. above),
- two-step correction impression (2-step/2-phase) with R-SI-LINE® PUTTY or R-SI-LINE® PUTTY FS,
- one-phase (monophase) impressions for crowns, bridges, inlays, partial and complete dentures and
- relining impressions, functional impressions and impressions for prosthetic restorations.

### Performance features

The performance features of the product meet the requirements of the intended use.

### Contraindications

Not known.

### Patient target group

Persons who are treated during a dental procedure.

### Intended users

This medical device should only be used by a professionally trained dental practitioner.

### Application

#### 1. Preparing the auto-mixing cartridge/R-SI-LINE MIXING PISTOL N

Introduce the auto-mixing cartridge into the MIXING PISTOL N (REF MIP1081) and lock it. Before first use screw off the original cap and dispose it. When screwing on a new mixing cannula the cartridge should be tested if both pastes will be extruded equally. Then screw on the mixing cannula R-SI-LINE MIXING CANNULA. The mixing pistol is ready for application.

For direct fine application of the R-SI-LINE® MEDIUM SH (- FS) onto the prepared areas additionally the R-SI-LINE® INTRA-ORAL-TIPS L should be used.

For fine applications instead of R-SI-LINE® MIXING CANNULA LN (REF MIK1099) + INTRA-ORAL-TIP L (REF IOT2001) a NEEDLE MIXING CANNULA MLNN (REF MIK3099N) with an integrated, curved application needle is available.

#### 2 a. One-step sandwich-technique

The impression tray will be loaded with the base impression material.

Before every application a small amount of material (size of a pea) should be pressed out and should be disposed immediately.

Load the auto-mixed R-SI-LINE® MEDIUM SH (- FS) partially/voluminous onto the teeth or the prepared areas and partially/voluminous onto the base impression material in the tray which is coated with R-SI-LINE® SILICONE ADHESIVE (REF SIA1070).

The process should begin in such a way that the application around the preparation and the filling of the tray are completed at the same time.

Insert the tray with R-SI-LINE® MEDIUM SH into the mouth not later than 1 min 30 s (the fast set component R-SI-LINE® MEDIUM SH FS not later than 60 s) after extruding (working time at 74 °C).

Maintain the position for R-SI-LINE® MEDIUM SH at least 3 min (the fast set R-SI-LINE® MEDIUM SH FS at least 2 min) in the mouth.

Note: The filled MIXING CANNULA functions as a seal until the next application. Closing the cartridge by turning on the original cap between use could cause cross contamination of base - and catalyst paste and thus premature setting.

#### 2 b. Two-step correction impression technique

First cut drain channels into the pre-impression.

Before every application a small amount of material (size of a pea) should be pressed out and should be disposed immediately.

Load the auto-mixed R-SI-LINE® MEDIUM SH (- FS) directly onto the teeth or the prepared areas and partially/voluminous onto the base impression material (R-SI-LINE® PUTTY, - PUTTY FS, - PUTTY MATIC/ - HEAVY MATIC, - HEAVY KINGSIZE) in the tray which is coated with R-SI-LINE® SILICONE ADHESIVE (REF SIA1070).

Insert the tray with R-SI-LINE® MEDIUM SH into the mouth not later than 2 min (the fast set component R-SI-LINE® MEDIUM SH FS not later than 1 min 30 s) after extruding (working time at 74 °C).

Maintain the position for R-SI-LINE® MEDIUM SH at least 3 min (the fast set R-SI-LINE® MEDIUM SH FS at least 2 min) in the mouth.

Note: The filled MIXING CANNULA functions as a seal until the next application. Closing the cartridge by turning on the original cap between use could cause cross contamination of base - and catalyst paste and thus premature setting.

#### 3. Casting

The impression can be poured 30 min after removal from the mouth. There are no other time limitations.

#### 4. Electroplating

All impressions can be copper- or silverplated.

#### 5. Disinfection

The impression can be disinfected with disinfectant solutions especially for materials based on addition-curing silicones. Please follow the manufacturer's instructions for use.

#### 6. Storage

Do not store above 25 °C (77 °F)! Avoid exposure to sunlight. Do not use after expiry date.

#### Additional notes

Do not use R-SI-LINE® MEDIUM SH (- FS) with standard condensation reaction type silicone impression materials.

Do not use latex or sulfur-containing polymeric gloves because they may interfere with the setting reaction of the impression material.

When using the sandwich-technique (1-step/2-phase), make sure that both materials have the same setting time in mouth. The working time is determined by the material with the shorter working time.

Keep out of the reach of children!

#### Disposal

Disposal of the product according to local authority regulations.

#### Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

R-SI-LINE is a registered trademark of R-dental Dentalerzeugnisse GmbH.

#### Warranty

R-dental Dentalmaterialien GmbH warrants this product will be free from defects in material and manufacture. R-dental Dentalmaterialien GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. The user is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and R-dental Dentalmaterialien GmbH's sole obligation shall be repair or replacement of the product.

#### Limitation of Liability

Except where prohibited by law, R-dental Dentalmaterialien GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

© 2024 R-dental Dentalmaterialien GmbH. All rights reserved.